


MiniMed™ Mio™

Infusion set • Cathéter • Infusionsset • Equipo de infusión • Infusieset • Set di infusione • Slangesett • Infusiosetti • Infusionsset • Infusionsæt • Conjunto de infusão • Infuzijski set • Infuzní set • Infúziós szerelék • Zestaw infuzyjny • Infúzyon seti • Infúzna súprava • Σετ έγχυσης • Инфузионный набор • 输注管路 • أداة التثريب • טט עירוי





Medtronic


Front • Avant • Vorderseite • Cuerpo de frente •
 Voorkant • Davanti • Foran • Edessä • Fram • Forside •
 Frente • Spredaj • Čelní strana • előlnézet • Przód • Ön •
 Vpredu • Εμπρός • Живот • 正面 • الأمام • 正門

Back • Arrière • Rückseite • Cuerpo de espalda •
 Achterkant • Dietro • Bak • Takana • Bak • Bagside •
 Costas • Zadaj • Zadní strana • hátulnézet • Tyl • Arka •
 Vzadu • Плов • Спина • 背面 • الخلف • 背


Click • Clic • Klick • Hacer clic • Vastklikken • Clic •
 Kliikk • Napsahdus • Klick • Klik • Clique • Klik •
 Zacvaknout • kattanas • Kliknięcie • Tik • Cvaknutie •
 Клік • Щелчок • 咔嚓一声 • طقة • קליק


4


5


6


7


12


13

14



15

Click • Clic • Klick • Hacer clic • Vastklikken • Clic •
 Klikk • Napsahdus • Klick • Klik • Clique • Klik •
 Zacvaknout • kattanás • Kliknięcie • Tik • Cvaknutie •
 Κλικ • Щелчок • 咔嚓一声 • قطة • קליק




DISCONNECTING • DÉCONNEXION • A
 TRENNEN • DESCONEXIÓN •
 LOSKOPPELEN • SCOLLEGAMENTO •
 KOBLE FRA • IRROTTAMINEN •
 BORTKOPPLING • FRAKOBLING •
 DESLIGAR • ODKLÓP • ODPOJENÍ •
 LEVÁLASZTÁS • ODŁĄCZANIE •
 BAĞLANTIYI KESME • ODPOJENIE •
 ΑΠΟΣΥΝΔΕΣΗ • ОТСОЕДИНЕНИЕ •
 断开连接 • الفصل • ניתוק


Click • Clic • Klick • Hacer clic • Vastklikken • Clic • Klick •
 Napsahdus • Klick • Klik • Clique • Klik • Zacvaknout • kattánás •
 Kliknięcie • Tik • Cvaknutie • Κλικ • Щелчок • 咔嚓一声 • طقة • קליק

RECONNECTING • RECONNEXION •
 WIEDERANSCHLIESSEN •
 RECONEXIÓN • AANSLUITEN •
 RICOLLEGAMENTO • KOBLE TIL
 IGJEN • KYTKEMINEN TAKAISIN •
 ÅTERANSLUTNING •
 GENTILKOBLING • VOLTAR A LIGAR •
 PONOVNÍ PRIKLOP • OPĚTOVNÉ
 PŘIPOJENÍ •
 ÚJRACSATLAKOZTATÁS • PONOWNE
 PODŁĄCZANIE • YENİDEN
 BAĞLAMA • OPĀTOVNÉ PRIPOJENIE •
 ΕΠΑΝΑΣΥΝΔΕΣΗ • ПОВТОРНОЕ
 ПОДСОЕДИНЕНИЕ • 重新连接 •
 חיבור מחדש • إعادة التوصيل

Click • Clic • Klick • Hacer clic •
 Vastklikken • Clic • Klick • Napsahdus •
 Klick • Klik • Clique • Klik • Zacvaknout •
 kattánás • Kliknięcie • Tik • Cvaknutie •
 Κλικ • Щелчок • 咔嚓一声 • طقة • קליק

MiniMed™ Mio™

EN See page 8

FR Consultez la page 12

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ES Consultar la página 20

NL Zie blz. 24

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NO Se side 32

FI Katso sivu 36

SV Se sidan 40

DA Se side 44

PT Consultar a página 48

SL Glejte stran 52

CS Viz strana 56

HU Lásd: 60. oldal

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TR Bkz. sayfa 68

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Indications for use

The MiniMed Mio infusion set is indicated for the subcutaneous infusion of insulin from an infusion pump. The infusion set is indicated for single use.

Intended use

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulins approved for subcutaneous infusion.

Description

The MiniMed Mio is an infusion device which combines an infusion set (catheter) with an aid for insertion (serter). As the infusion set is an integrated part of the serter, the infusion set is delivered assembled and ready for use as a single, sterile unit.

Contraindications

The infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous (IV) infusion. Do not use the infusion set for infusion of blood or blood products.

Warnings

- Read all instructions carefully before using the infusion set. Failure to follow instructions may result in pain or injury.
- The infusion set is sterile and non-pyrogenic unless the package has been opened or damaged. Do not use the infusion set if the package has been opened or damaged.

- For Medtronic pump users, do not use the infusion set if the tubing connector needle has been damaged.
- Ensure sterility by checking that the sterile paper and tamper-proof seal are not damaged.
- If using this infusion set for the first time, do the first set-up in the presence of a healthcare professional.
- Inaccurate medication delivery, site irritation or infection may result from improper insertion and maintenance of the infusion site.
- Do not re-insert the introducer needle into the infusion set. This can cause tearing or damage to the soft cannula and can lead to unpredictable medication flow.
- Reuse of the infusion set can cause site irritation or infection and can lead to inaccurate medication delivery.
- Always remove the needle guard before inserting the infusion set.
- Before insertion, clean the insertion site with a disinfectant as indicated by a healthcare professional.
- For the MiniMed tubing connector: If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.

Precautions

- Never try to fill or free clogged tubing while the infusion set is inserted. This can result in unpredictable medication flow.
- When unwinding the tubing, release it with caution, because a hard pull on the tubing can result in damage to the infusion set or introducer needle. Ensure that the infusion set is properly in place when the tubing is fully released.
- Do not leave air in the infusion set. Make sure to fill the tubing completely. Consult the pump manufacturer's instructions on how to fill the tubing.
- The infusion set is a single use item, which must be disposed of after use. Do not clean or re-sterilize.
- Do not put disinfectants, perfumes, or deodorants on the infusion set as these may impact the integrity of the infusion set.
- Replace the infusion set if the adhesive tape becomes loose, or if the soft cannula is displaced from its original site.
- Check the infusion site frequently to ensure that the soft cannula remains firmly in place. If the infusion set is not in place, replace it with a new one at a new site. Since the cannula is soft, it will not cause any pain if it slips out, which may happen without notice. The soft cannula must always be completely inserted to receive the full amount of medication.
- Replace the infusion set every two to three days, or per a healthcare professional's instructions.

- If the infusion site becomes irritated or infected, replace the set, and use a new site until the first site has healed.
- Use aseptic techniques when temporarily disconnecting the infusion set. Consult a healthcare professional for how to compensate for missed medication while disconnected.
- Carefully monitor blood glucose levels when disconnected, and after reconnecting the infusion set.
- Never point a loaded insertion device towards any body part, where insertion is not desired.
- Check blood glucose levels one to three hours after inserting the infusion set. Measure blood glucose on a regular basis. Discuss this with a healthcare professional.
- Do not change the infusion set just before bedtime unless blood glucose can be checked one to three hours after insertion.
- If blood glucose becomes inexplicably high, or an occlusion alarm occurs, check for clogs or leaks. If in doubt, change the infusion set, because the soft cannula may be dislodged or partially clogged. Discuss a plan to lower blood glucose levels with a healthcare professional. Test blood glucose to ensure the problem has been corrected.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and your national authority.

Storage and disposal

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Put the lid back on the insertion device and dispose in a proper sharps container according to local laws, to avoid risk of needle stick injury.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

Warranty


- For product warranty information, please contact a local Medtronic support representative or visit: www.medtronicdiabetes.com/warranty.

© 2020 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (TM*) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.

MiniMed™ Mio™

Preparation and insertion

This product is an infusion device which combines an infusion set (catheter) with an aid for insertion (serter). This single-use, sterile product is delivered assembled and ready for use. Ensure sterility by checking that the sterile paper and tamper-proof seal are not damaged. Carefully examine the instructions for use before inserting this product. Refer to pump and reservoir user guides for installation and filling. Please see the corresponding illustrations on pages 2 and 3.

Contents

- A. Lid
- B. Needle guard
- C. Soft cannula
- D. Adhesive tape
- E. Cannula housing
- F. Site connector
- G. Tubing
- H. Introducer needle
- I. Insertion device
- J. Sterile paper

Important note!

- K. Three raised dots – press here when lifting the lid off (see Figure 6)
- L. Lined indentations – press here when loading (see Figure 10)
- M. Round Indentations – press both sides when releasing (see Figure 12)
- N. Disconnect cover (see Figure B)

Instructions

1. Wash your hands.
2. Suggested infusion sites. Choose infusion sites recommended by a healthcare professional. Do not use an area right next to your previous insertion site.
3. Clean the insertion area with a disinfectant as directed by your healthcare professional. Be sure the area is air dry before inserting the infusion set.
4. Pull the red tape to remove the protective seal.
5. Remove the sterile paper.
6. Hold the bottom of the insertion device with one hand, while firmly pressing the three raised dots on each side of the lid with the other hand, and lift up the lid.
CAUTION: Be careful not to bend or touch the introducer needle when preparing theserter and infusion set.
7. **Unwinding the tubing:** Carefully twist the beginning of the tubing out of the slot. Now unwind the tubing from the insertion device by pulling it carefully upwards.
Caution: Do not touch the introducer needle when unwinding the tubing. Be careful not to pull too hard when the tubing is fully released. You could accidentally remove the infusion set from the introducer needle.
IMPORTANT NOTE! If the infusion set is not securely placed in the insertion device with the needle pointing straight ahead prior to insertion, pain or minor injury may occur. To avoid this, please ensure that the following steps are taken

gently so the infusion set is not loosened / removed by accident.

8. Connect the tubing to a filled reservoir. Fill the infusion set. Ensure that there are no air bubbles in the tubing. Consult the pump manufacturer's instructions for use.

Caution: When filling the infusion set, hold the infusion set with the needle pointing down to ensure that medication does not come in contact with the adhesive backing paper.

9. Pull up gently to remove the backing paper.
10. **To prepare the serter:** Place fingers on lined indentations. Pull up the spring until you hear a "click".
11. Carefully remove the needle guard by gently twisting and then pulling it off. Make sure that the soft cannula does not exceed the introducer needle.

Caution: Before positioning the infusion set over an insertion site, place the tubing in the slot on the side of the serter to ensure that the tubing is not caught underneath the infusion set during insertion.
12. **To insert the infusion set:** Position the infusion set over the insertion site. Simultaneously press the round indentations on both sides of the insertion device to insert the infusion set.
13. Push gently on the spring to secure the adhesive onto the skin.
14. Remove the serter and introducer needle by grasping the center of the serter and pulling it gently back. Gently massage the adhesive tape

onto the skin. Fill the empty space in the cannula with the following amount of U-100 insulin:

U-100 insulin:

6 mm: 0.3 units (0.003 ml)

9 mm: 0.5 units (0.005 ml)

15. Push the lid back in place until you hear a "click."

To disconnect the infusion set

The infusion set allows you to temporarily disconnect from your pump without having to change the infusion set.

- A. Hold the adhesive tape in place by placing a finger just in front of the cannula housing, and gently squeeze the sides of the site connector. Pull the site connector straight out from the cannula housing.
- B. Insert the disconnect cover into the cannula housing until you hear a "click".

To reconnect the infusion set

- C. If necessary, fill the infusion set until medication dispenses from the needle.

Caution: When filling the infusion set, hold the needle pointing down.
- D. Remove the disconnect cover from the cannula housing by pressing the grips. Place a finger on the adhesive tape just in front of the cannula housing.
- E. Push the connector needle straight in until you hear a "click"

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Medtronic Ltd.
Medtronic Österreich GmbH